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OPERATING TECHNIQUE

“BASIC-NAIL”

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
.....

*Locked nail for proximal fractures
of the femur in Ti
uncannulated*

- *Modular system, standard, medium*
- *Reduced valgus angle*
- *Dynamic or static implant*


fig. 1

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• *Introduce the nail from the apex of the great trochanter through the common supra-trochanteric access and then, under X-ray control, prepare an entry for the diaphyseal endo-midullary component (nail) by the trochanteric ridge.*


fig. 2

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• *Having decided on the point of introduction of the nail.*

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• *Use the Aul (14-0040483) to perforate, checking both views with the image amplifier as you go.*

fig. 3


fig. 4

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• *Manually introduce the nail fixed to its inserter into the medullary cavity of the femur. Drive in the nail till the axis of the lag screw is positioned in the lower quadrant of the femur neck. The nail should advance without difficulty by simply pushing with your hand.*

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• *Should this insertion require a slight force, use the Impactor (14-0040495), which must be placed on top of the nail fixing bolt.*

▲ **Carry out an X-ray check.**

▲ **Avoid any kind of percussion against the nail inserter.**


fig. 6


fig. 7

fig. 8

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• *Insert the Lag Screw Guide Sleeve (14-0040430) in the nail inserter till it touches the bone. The lag screw guide sleeve must stay in place fixed to the nail inserter with its knob until just before the introduction of the distal screws; introduce the Proximal awl (14-0040440) and then penetrate the lateral cortical.*

▲ **Fix the sleeve to the nail inserter using its knob to stabilize the system.**

▲ **Don't tighten the knob excessively as this may deform the sleeve.**

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• Remove the awl and insert the Cannulated Proximal Probe (14-0040435) into the guide sleeve with the 3x410mm Guide Wire (14-0040445). The probe should be pushed in till it touches the bone.

Push in the wire through the lateral cortical till it reaches the subchondral zone of the femur head.

▲ Carry out X-ray checks (anteroposterior and axial).

▲ If the X-rays show the guide wire is bent, it must absolutely be removed and re-inserted correctly or even replaced with another one.


fig. 10

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• Having correctly positioned the guide wire, measure the part of it remaining outside the probe by placing it in the slot of the Lag Screw Gauge (14-0030575) resting against the end of the probe itself.


fig. 11

BASIC-NAIL LOCKED NAIL - OPERATING TECHNIQUE

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• *Extract the probe from the guide sleeve being careful that the guide wire too doesn't move back. Introduce the Lag Screw Drill (14-0040484).*

⚠ **Drill until coming into contact with the sleeve, thus obtaining a housing for the lag screw.**


fig. 12

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• *Having extracted the lag screw drill keeping the guide wire well in place, put the screw in through the lag screw guide sleeve. To do this, use the Lag Screwdriver (14-0040485).*

Screw in till the symbol marked on the stem of the screwdriver reaches the edge of the sleeve.

Twist the ring clockwise as far as the stop and proceed with compression where necessary.


fig. 13

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• *Insert the Distal Guide Sleeve (14-0030585) in the chosen hole in the nail inserter till it touches the cortical. "Round" for the round hole, "Stat" (static) for the proximal part of the slot and "Din" (dynamic) for distal part of the slot. If the 250mm nail is used, choose the holes provided for that nail.*

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• *Put into the sleeve the 3,5 mm Distal Guide Sleeve (14-0040460) with the Distal Drill Bit (14-0040480) inside.*

⚠ **Perforate the two corticals advancing slowly.**

⚠ **Carry out X-ray checks.**

BASIC-NAIL LOCKED NAIL - OPERATING TECHNIQUE

fig. 14

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- The Basic-Nail locked femoral nail (Standard or 250 mm version) is an alternative to open synthesis techniques. It provides reduction and immobilization of the fracture by a closed technique which limits the risk of infection and blood loss. As is normal with all orthopaedic implant systems, there is a risk of failure and this can vary from patient to patient.

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Indications and contraindications

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- Choosing this type of implant and the size to use is a matter of general traumatological principles and of the surgeon's experience. Generally speaking, the Basic-Nail is used for:

- 1) basicervical fractures;
- 2) pertrochanteric fractures;
- 3) intertrochanteric fractures;
- 4) sub-trochanteric fractures.

The implant should not be used when there is:

- inadequate quantity and quality of bone tissue;
- infection in progress;
- excessive geometrical anomaly in the shape of the femur.

As with all osteosynthesis systems and in spite of correct sizing, there is a risk of fatigue-induced breakage in the event of delayed-union, non-union or movement following removal of the distal screws.

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- Evo-Nail and Evo-Nail 250 mm locked femoral nails and all their components are made of titanium alloy, whose composition and mechanical characteristics conform to national and international standards (ISO 5832/3).

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- All components of the implant must be handled and stored with extreme care, as any alteration (deformation, scratches, and so on) could put the implant at risk. In the operating theatre it is likewise necessary to make sure the nail isn't impaired by the instruments. For reasons of compatibility between materials and geometrical form, only GRUPPO BIOIMPIANTI lag screws, blades and distal screws should be used. As with implant materials, great care must also be taken with the instruments, which must be stored in their racks to prevent accidental impact and surface damage.

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- Each nail is supplied appropriately packaged, marked and sterile by β rays.

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GRUPPO BIOIMPIANTI S.r.l.