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*"Esoblock®"
Reabsorbable diaphyseal plug*

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The ESOBLOCK® biodegradable diaphyseal plug enables the medullary cavity to be sealed off before acrylic cement is introduced.

It stops the cement moving inside in the diaphysis, facilitates the pressurization of the cement and guarantees excellent protection of the intramedullary tissue.

Being perfectly biocompatible, it is re-absorbed in 1 to 3 days at 37° and withstands the momentary increase in temperature due to polymerization of the cement.

MODULAR DESIGN AND COMPOSITION

Each plug comes in double, beta-ray sterilized packaging, in the following sizes: 8, 10, 12, 14, 16, 18 and 20, the last two available on special request. Each size matches exactly the diameter of the medullary canal. For example, size 10 provides optimum obturation for a 10 mm canal diameter.

The plug is made of: porcine gelatine, glycerol, water, methylparahydroxybenzoate.

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CHARACTERISTICS

- ESOBLOCK® diaphyseal plugs are totally re-absorbable, facilitate subsequent prosthetic revision and, unlike synthetic plugs, they minimize loss of bone stock.
- Their gelatine based chemical composition gives them elastic properties that ensure optimum filling of the medullary cavity as well as improved resistance to the pressurization of the cement. The plug takes the anatomical shape of the canal. Synthetic plugs, on the other hand, do not allow such deformation as they are more rigid and cannot adhere perfectly to the canal or guarantee a seal against eventual movement of the cement.
- None of the "debris" usually found when using synthetic plugs.
- Easy to implant.

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- 1 Fit the trial plug on the insert device and proceed with the measurement (fig. A).
- 2 Position the cursor at the desired distance (distance between the top edge of the bone and the distal end of the prosthesis to implant); add 1 cm to this distance (fig. B).
- 3 Determine the plug size by measuring the diameter of the femoral canal with the trial plug.
- 4 Remove the trial plug and mount the femur plug holder.
- 5 Put the femur plug in the plug holder and proceed to insert it in the medullary canal (fig. C).
- 6 Pull out the insert device and introduce the acrylic cement.